

Product Datasheet

Human Calponin 1 ELISA Kit (orb1213941)

Description	Human Calponin 1, CNN1 ELISA Kit
Reactivity	Human
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	CNN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	33,170
Uniprot ID	P51911
Solubility (25°C)	16ng/ml
Sensitivity	0.02ng/ml
Entrez	1264
Expiration Date	Please enquire.