

Product Datasheet

Human Mothers Against Decapentaplegic Homolog 3 ELISA Kit (orb1213957)

Description	Human Mothers Against Decapentaplegic Homolog 3, SMAD3 ELISA Kit
Reactivity	Human
Range	0.1-35ng/ml
Tested Applications	ELISA
Target	SMAD3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	48,081
Uniprot ID	P84022
Solubility (25°C)	40ng/ml
Sensitivity	0.047ng/ml
Entrez	4088
Expiration Date	Please enquire.