

Product Datasheet

Human Mothers Against Decapentaplegic Homolog 2 ELISA Kit (orb1213958)

Description	Human Mothers Against Decapentaplegic Homolog 2, SMAD2 ELISA Kit
Reactivity	Human
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	SMAD2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	52,306
Uniprot ID	Q15796
Solubility (25°C)	12.8ng/ml
Sensitivity	0.019ng/ml
Entrez	4087
Expiration Date	Please enquire.