

Product Datasheet

Human Ubiquilin 2 ELISA Kit (orb1214058)

Description	Human Ubiquilin 2, UBQLN2 ELISA Kit
Reactivity	Human
Range	50-10000ng/L
Tested Applications	ELISA
Target	UBQLN2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	65,696
Uniprot ID	Q9UHD9
Solubility (25°C)	12800ng/L
Sensitivity	25.26ng/L
Entrez	29978
Expiration Date	Please enquire.