

Product Datasheet

Human Metalloreductase Steap2 ELISA Kit (orb1214069)

| | |
|----------------------------|--|
| Description | Human Metalloreductase Steap2, STEAP2 ELISA Kit |
| Reactivity | Human |
| Range | 20-7000ng/L |
| Tested Applications | ELISA |
| Target | STEAP2 |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 56,056 |
| Uniprot ID | Q8NFT2 |
| Solubility (25°C) | 8000ng/L |
| Sensitivity | 9.96ng/L |
| Entrez | 261729 |
| Expiration Date | Please enquire. |