

Product Datasheet

Human Proliferation Associated Protein 2G4 ELISA Kit (orb1214080)

Description	Human Proliferation Associated Protein 2G4, PA2G4 ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	PA2G4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	43,787
Uniprot ID	Q9UQ80
Solubility (25°C)	4000ng/L
Sensitivity	8.77ng/L
Entrez	5036
Expiration Date	Please enquire.