

Product Datasheet

Human Ghrelin-O-Acyltransferase ELISA Kit (orb1214120)

Description	Human Ghrelin-O-Acyltransferase, GOAT ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	MBOAT4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	49,716
Uniprot ID	Q96T53
Solubility (25°C)	64ng/ml
Sensitivity	0.12ng/ml
Entrez	619373
Expiration Date	Please enquire.