

Product Datasheet

Human Growth Differentiation Factor 5 ELISA Kit (orb1214168)

Description	Human Growth Differentiation Factor 5, GDF5 ELISA Kit
Reactivity	Human
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	GDF5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	55,411
Uniprot ID	P43026
Solubility (25°C)	24ng/ml
Sensitivity	0.026ng/ml
Entrez	8200
Expiration Date	Please enquire.