

Product Datasheet

Human Glycyl Trna Synthetase ELISA Kit (orb1214180)

Description	Human Glycyl Trna Synthetase, GARS ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	SLC27A5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	75,385
Uniprot ID	Q9Y2P5
Solubility (25°C)	4800ng/L
Sensitivity	10.03ng/L
Entrez	10998
Expiration Date	Please enquire.