

Product Datasheet

Human Leukocyte Elastase Inhibitor ELISA Kit (orb1214182)

Description	Human Leukocyte Elastase Inhibitor, LEI ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	SERPINB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,742
Uniprot ID	P30740
Solubility (25°C)	48ng/ml
Sensitivity	0.057ng/ml
Entrez	1992
Expiration Date	Please enquire.