

Product Datasheet

Human Glial Fibrillary Acidic Protein ELISA Kit (orb1214207)

Description	Human Glial Fibrillary Acidic Protein, GFAP ELISA Kit
Reactivity	Human
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	GFAP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	49,880
Uniprot ID	P14136
Solubility (25°C)	16ng/ml
Sensitivity	0.026ng/ml
Entrez	2670
Expiration Date	Please enquire.