

Product Datasheet

Human Chitinase-3-like Protein 1 ELISA Kit (orb1214215)

Description	Human Chitinase-3-like Protein 1, CHI3L1 ELISA Kit
Reactivity	Human
Range	1-400ng/ml
Tested Applications	ELISA
Target	CHI3L1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,625
Uniprot ID	P36222
Solubility (25°C)	480ng/ml
Sensitivity	0.61ng/ml
Entrez	1116
Expiration Date	Please enquire.