

Product Datasheet

Human Lipocalin 1 ELISA Kit (orb1214217)

Description	Human Lipocalin 1, LCN1 ELISA Kit
Reactivity	Human
Range	0.5-200ng/mL
Tested Applications	ELISA
Target	LCN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	19,250
Uniprot ID	P31025
Solubility (25°C)	240ng/mL
Sensitivity	0.35ng/mL
Entrez	3933
Expiration Date	Please enquire.