

Product Datasheet

Human Bone Morphogenetic Protein 4 ELISA Kit (orb1214228)

Description	Human Bone Morphogenetic Protein 4, BMP-4 ELISA Kit
Reactivity	Human
Range	1-400ng/ml
Tested Applications	ELISA
Target	BMP4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,555
Uniprot ID	P12644
Solubility (25°C)	480ng/ml
Sensitivity	0.57ng/ml
Entrez	652
Expiration Date	Please enquire.