

Product Datasheet

Human Beta-defensin 1 ELISA Kit (orb1214238)

Description	Human Beta-defensin 1, DEFB1 ELISA Kit
Reactivity	Human
Range	5-2000pg/ml
Tested Applications	ELISA
Target	DEFB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	7,420
Uniprot ID	P60022
Solubility (25°C)	2400pg/ml
Sensitivity	2.31pg/ml
Entrez	1672
Expiration Date	Please enquire.