

Product Datasheet

Human Intercellular Adhesion Molecule 4 ELISA Kit (orb1214239)

Description	Human Intercellular Adhesion Molecule 4, ICAM-4 ELISA Kit
Reactivity	Human
Range	2-600ng/ml
Tested Applications	ELISA
Target	ICAM4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,265
Uniprot ID	Q14773
Solubility (25°C)	800ng/ml
Sensitivity	1.12ng/ml
Entrez	3386
Expiration Date	Please enquire.