

Product Datasheet

Human High Mobility Group Protein B1 ELISA Kit (orb1214273)

Description	Human High Mobility Group Protein B1, HMGB-1 ELISA Kit
Reactivity	Human
Range	0.5-150ng/ml
Tested Applications	ELISA
Target	HMGB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,894
Uniprot ID	P09429
Solubility (25°C)	160ng/ml
Sensitivity	0.24ng/ml
Entrez	3146
Expiration Date	Please enquire.