

Product Datasheet

Human Substance P Receptor ELISA Kit (orb1214294)

Description	Human Substance P Receptor, SP-R ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	TACR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,251
Uniprot ID	P25103
Solubility (25°C)	64ng/ml
Sensitivity	0.12ng/ml
Entrez	6869
Expiration Date	Please enquire.