

Product Datasheet

Human Cholecystokinin ELISA Kit (orb1214319)

Description	Human Cholecystokinin, CCK ELISA Kit
Reactivity	Human
Range	2-100ng/L
Tested Applications	ELISA
Target	CCK
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	12,669
Uniprot ID	P06307
Solubility (25°C)	128ng/L
Sensitivity	1.32ng/L
Entrez	885
Expiration Date	Please enquire.