

Product Datasheet

Human Motilin ELISA Kit (orb1214320)

Description	Human Motilin, MTL ELISA Kit
Reactivity	Human
Range	5-1500ng/L
Tested Applications	ELISA
Target	MYOT
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	55,395
Uniprot ID	Q9UBF9
Solubility (25°C)	1600ng/L
Sensitivity	2.51ng/L
Entrez	9499
Expiration Date	Please enquire.