

Product Datasheet

Human Retinol Binding Protein 4 ELISA Kit (orb1214339)

Description	Human Retinol Binding Protein 4, RBP-4 ELISA Kit
Reactivity	Human
Range	1-400ng/ml
Tested Applications	ELISA
Target	RBP4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,010
Uniprot ID	P02753
Solubility (25°C)	480ng/ml
Sensitivity	0.51ng/ml
Entrez	5950
Expiration Date	Please enquire.