

Product Datasheet

Human Platelet Basic Protein ELISA Kit (orb1214351)

Description	Human Platelet Basic Protein, PBP ELISA Kit
Reactivity	Human
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	PPBP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	13,894
Uniprot ID	P02775
Solubility (25°C)	48ng/mL
Sensitivity	0.043ng/mL
Entrez	5473
Expiration Date	Please enquire.