

Product Datasheet

Human Alpha1 Antichymotrypsin ELISA Kit (orb1214365)

Description	Human Alpha1 Antichymotrypsin, AACT ELISA Kit
Reactivity	Human
Range	10-3500ug/ml
Tested Applications	ELISA
Target	SERPINA3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	47,651
Uniprot ID	P01011
Solubility (25°C)	4000µg/ml
Sensitivity	5.17µg/ml
Entrez	12
Expiration Date	Please enquire.