

Product Datasheet

Human Hepcidin ELISA Kit (orb1214375)

Description	Human Hepcidin, HEPC ELISA Kit
Reactivity	Human
Range	10-4000pg/ml
Tested Applications	ELISA
Target	HEPC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	9,408
Uniprot ID	P81172
Solubility (25°C)	4800pg/ml
Sensitivity	5.12pg/ml
Entrez	57817
Expiration Date	Please enquire.