

Product Datasheet

Human 5-Nucleotidase ELISA Kit (orb1214408)

Description	Human 5-Nucleotidase, 5-NT ELISA Kit
Reactivity	Human
Range	0.5-150ng/mL
Tested Applications	ELISA
Target	NT5E
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	63,368
Uniprot ID	P21589
Solubility (25°C)	160ng/mL
Sensitivity	0.34ng/mL
Entrez	4907
Expiration Date	Please enquire.