

Product Datasheet

Human Indoleamine 2 ELISA Kit (orb1214412)

Description	Human Indoleamine 2, 3-dioxygenase, IDO ELISA Kit
Reactivity	Human
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	IDO1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	45,326
Uniprot ID	P14902
Solubility (25°C)	96ng/ml
Sensitivity	0.16ng/ml
Entrez	3620
Expiration Date	Please enquire.