

Product Datasheet

Human Dynamin 1 ELISA Kit (orb1214438)

Description	Human Dynamin 1, DNM1 ELISA Kit
Reactivity	Human
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	DNM1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	97,408
Uniprot ID	Q05193
Solubility (25°C)	48ng/mL
Sensitivity	0.044ng/mL
Entrez	1759
Expiration Date	Please enquire.