

Product Datasheet

Human Mothers Against Decapentaplegic Homolog 1 ELISA Kit (orb1214458)

Description	Human Mothers Against Decapentaplegic Homolog 1, SMAD1 ELISA Kit
Reactivity	Human
Range	0.5-150ng/ml
Tested Applications	ELISA
Target	SMAD1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	52,260
Uniprot ID	Q15797
Solubility (25°C)	160ng/ml
Sensitivity	0.25ng/ml
Entrez	4086
Expiration Date	Please enquire.