

Product Datasheet

Human Aquaporin 4 ELISA Kit (orb1214460)

Description	Human Aquaporin 4, AQP-4 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	AQP4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	34,830
Uniprot ID	P55087
Solubility (25°C)	24ng/ml
Sensitivity	0.022ng/ml
Entrez	361
Expiration Date	Please enquire.