

Product Datasheet

Human Insulin-like Growth Factors Binding Protein 4 ELISA Kit (orb1214473)

Description	Human Insulin-like Growth Factors Binding Protein 4, IGFBP4 ELISA Kit
Reactivity	Human
Range	30-9000ng/mL
Tested Applications	ELISA
Target	IGFBP4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,934
Uniprot ID	P22692
Solubility (25°C)	9600ng/mL
Sensitivity	14.54ng/mL
Entrez	3487
Expiration Date	Please enquire.