

Product Datasheet

Human Insulin-like Growth Factors Binding Protein 3 ELISA Kit (orb1214476)

Description	Human Insulin-like Growth Factors Binding Protein 3, IGFBP-3 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	IGFBP3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	31,674
Uniprot ID	P17936
Solubility (25°C)	24ng/ml
Sensitivity	0.026ng/ml
Entrez	3486
Expiration Date	Please enquire.