

## Product Datasheet

# Human Complement Factor B ELISA Kit (orb1214478)

<b>Description</b>	Human Complement Factor B, CFB ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	8-2000ug/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	CFB
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	85,533
<b>Uniprot ID</b>	<b>P00751</b>
<b>Solubility (25°C)</b>	2400µg/ml
<b>Sensitivity</b>	4.19µg/ml
<b>Entrez</b>	<b>629</b>
<b>Expiration Date</b>	Please enquire.