

## Product Datasheet

### Human Toll-like Receptor 1 ELISA Kit (orb1214483)

|                            |                                            |
|----------------------------|--------------------------------------------|
| <b>Description</b>         | Human Toll-like Receptor 1, TLR1 ELISA Kit |
| <b>Reactivity</b>          | Human                                      |
| <b>Range</b>               | 0.2-60ng/ml                                |
| <b>Tested Applications</b> | ELISA                                      |
| <b>Target</b>              | TLR1                                       |
| <b>Storage</b>             | 2-8°C                                      |
| <b>Note</b>                | For research use only                      |
| <b>Application notes</b>   | Sandwich                                   |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates     |
| <b>MW</b>                  | 90,291                                     |
| <b>Uniprot ID</b>          | <a href="#">Q15399</a>                     |
| <b>Solubility (25°C)</b>   | 80ng/ml                                    |
| <b>Sensitivity</b>         | 0.09ng/ml                                  |
| <b>Entrez</b>              | <a href="#">7096</a>                       |
| <b>Expiration Date</b>     | Please enquire.                            |