

Product Datasheet

Human Toll-like Receptor 4 ELISA Kit (orb1214490)

Description	Human Toll-like Receptor 4, TLR4 ELISA Kit
Reactivity	Human
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	TLR4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	95,680
Uniprot ID	O00206
Solubility (25°C)	16ng/ml
Sensitivity	0.027ng/ml
Entrez	7099
Expiration Date	Please enquire.