

Product Datasheet

Human Triggering Receptor Expresses On Myeloid Cells-1 ELISA Kit (orb1214506)

Description	Human Triggering Receptor Expresses On Myeloid Cells-1, TREM-1 ELISA Kit
Reactivity	Human
Range	7-2000ng/L
Tested Applications	ELISA
Target	TREM1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,387
Uniprot ID	Q9NP99
Solubility (25°C)	2400ng/L
Sensitivity	3.59ng/L
Entrez	54210
Expiration Date	Please enquire.