

Product Datasheet

Human Vasculoar Endothelial Cell Growth Factor Receptor 2 ELISA Kit (orb1214511)

Description	Human Vasculoar Endothelial Cell Growth Factor Receptor 2, VEGFR-2 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	KDR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	151,527
Uniprot ID	P35968
Solubility (25°C)	48ng/ml
Sensitivity	0.05ng/ml
Entrez	3791
Expiration Date	Please enquire.