

Product Datasheet

Human Leptin Receptor ELISA Kit (orb1214515)

Description	Human Leptin Receptor, LEPR ELISA Kit
Reactivity	Human
Range	2-600ng/ml
Tested Applications	ELISA
Target	LEPR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	132,494
Uniprot ID	P48357
Solubility (25°C)	640ng/ml
Sensitivity	0.9ng/ml
Entrez	3953
Expiration Date	Please enquire.