

Product Datasheet

Human E-Cadherin ELISA Kit (orb1214523)

Description	Human E-Cadherin, E-CAD ELISA Kit
Reactivity	Human
Range	1-300ng/ml
Tested Applications	ELISA
Target	CDH1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	97,456
Uniprot ID	P12830
Solubility (25°C)	400ng/ml
Sensitivity	0.52ng/ml
Entrez	999
Expiration Date	Please enquire.