

Product Datasheet

Human Granulocyte Chemotactic Protein-2 ELISA Kit (orb1214525)

Description	Human Granulocyte Chemotactic Protein-2, GCP-2 ELISA Kit
Reactivity	Human
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	CXCL6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,897
Uniprot ID	P80162
Solubility (25°C)	12.8ng/ml
Sensitivity	0.022ng/ml
Entrez	6372
Expiration Date	Please enquire.