

Product Datasheet

Human Interleukin 1B ELISA Kit (orb1214535)

| | |
|----------------------------|---|
| Description | Human Interleukin 1B, IL-1B ELISA Kit |
| Reactivity | Human |
| Range | 20-6000pg/mL |
| Tested Applications | ELISA |
| Target | IL1B |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 30,748 |
| Uniprot ID | P01584 |
| Solubility (25°C) | 8000pg/mL |
| Sensitivity | 10.07pg/mL |
| Entrez | 3553 |
| Expiration Date | Please enquire. |