

Product Datasheet

Human Monocyte Chemotactic Protein 3 ELISA Kit (orb1214546)

Description	Human Monocyte Chemotactic Protein 3, MCP-3 ELISA Kit
Reactivity	Human
Range	5-1000ng/L
Tested Applications	ELISA
Target	CCL7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,200
Uniprot ID	P80098
Solubility (25°C)	1280ng/L
Sensitivity	3.07ng/L
Entrez	6354
Expiration Date	Please enquire.