

## Product Datasheet

# Human Macrophage Inflammatory Protein-3A ELISA Kit (orb1214568)

<b>Description</b>	Human Macrophage Inflammatory Protein-3A, MIP-3A ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	10-2000ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	CCL20
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	10,762
<b>Uniprot ID</b>	<b>P78556</b>
<b>Solubility (25°C)</b>	2400ng/L
<b>Sensitivity</b>	4.23ng/L
<b>Entrez</b>	<b>6364</b>
<b>Expiration Date</b>	Please enquire.