

## Product Datasheet

# Human Vascular Endothelial Cell Growth Factor C ELISA Kit (orb1214573)

<b>Description</b>	Human Vascular Endothelial Cell Growth Factor C, VEGF-C ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	20-6000ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	VEGFC
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	46,883
<b>Uniprot ID</b>	<a href="#">P49767</a>
<b>Solubility (25°C)</b>	6400ng/L
<b>Sensitivity</b>	1.18ng/L
<b>Entrez</b>	<a href="#">7424</a>
<b>Expiration Date</b>	Please enquire.