

Product Datasheet

Human Insulin ELISA Kit (orb1214600)

Description	Human Insulin, INS ELISA Kit
Reactivity	Human
Range	0.2-60mIU/L
Tested Applications	ELISA
Target	INS
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11981
Uniprot ID	P01308
Solubility (25°C)	80mIU/L
Sensitivity	0.11mIU/L
Entrez	3630
Expiration Date	Please enquire.