

Product Datasheet

Canine Coagulation Factor 8 ELISA Kit (orb1214648)

Description	Canine Coagulation Factor 8, F8 ELISA Kit
Reactivity	Canine
Range	1-400ng/ml
Tested Applications	ELISA
Target	F8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	265,831
Uniprot ID	O18806
Solubility (25°C)	480ng/ml
Sensitivity	0.59ng/ml
Entrez	403875
Expiration Date	Please enquire.