

## Product Datasheet

# Canine Metalloproteinase Inhibitor 1 ELISA Kit (orb1214656)

<b>Description</b>	Canine Metalloproteinase Inhibitor 1, TIMP1 ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	20-4500ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	TIMP1
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	22,852
<b>Uniprot ID</b>	<b>P81546</b>
<b>Solubility (25°C)</b>	4800ng/mL
<b>Sensitivity</b>	9.59ng/mL
<b>Entrez</b>	<b>403816</b>
<b>Expiration Date</b>	Please enquire.