

Product Datasheet

Canine Matrix Metalloproteinase 3 ELISA Kit (orb1214657)

Description	Canine Matrix Metalloproteinase 3, MMP-3 ELISA Kit
Reactivity	Canine
Range	0.5-200ng/mL
Tested Applications	ELISA
Target	MMP3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	53,633
Uniprot ID	Q6Y4Q5
Solubility (25°C)	240ng/mL
Sensitivity	0.27ng/mL
Entrez	403445
Expiration Date	Please enquire.