

## Product Datasheet

### Canine Interleukin 18 ELISA Kit (orb1214661)

<b>Description</b>	Canine Interleukin 18, IL-18 ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	6-1800ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	IL18
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	22,037
<b>Uniprot ID</b>	<a href="#">Q9XSR0</a>
<b>Solubility (25°C)</b>	1920ng/L
<b>Sensitivity</b>	3.12ng/L
<b>Entrez</b>	<a href="#">403796</a>
<b>Expiration Date</b>	Please enquire.