

Product Datasheet

Canine Coagulation Factor 5 ELISA Kit (orb1214677)

Description	Canine Coagulation Factor 5, F5 ELISA Kit
Reactivity	Canine
Range	0.5-200ng/ml
Tested Applications	ELISA
Target	F8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	265,831
Uniprot ID	O18806
Solubility (25°C)	240ng/ml
Sensitivity	0.26ng/ml
Entrez	403875
Expiration Date	Please enquire.