

Product Datasheet

Canine Erythropoietin ELISA Kit (orb1214680)

Description	Canine Erythropoietin, EPO ELISA Kit
Reactivity	Canine
Range	10-3000ng/L
Tested Applications	ELISA
Target	EPO
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,666
Uniprot ID	P33707
Solubility (25°C)	3200ng/L
Sensitivity	4.91ng/L
Entrez	404002
Expiration Date	Please enquire.